



PEPICO_Elettra_VG – Data Analysis in Python

Fabian Holzmeier (imec, Leuven, Belgium)

January 2022

Software Structure

Main program:






- `main.py` runs the code with the parameters provided in the input section

Classes:


- `VG_tof.py` reads in the raw files and creates lists with electron coordinates, ion stop times, start times, and assigns corresponding events
- `VG_tof_only.py` ion TOF spectra started by anode or random trigger
- `VG_figures.py` plots histograms, images, spectra
- `VG_Movie_Tof.py` produces electron energy selected TOF spectra and subtracts weighted TOF spectrum of randoms

Output: channel histograms

Channel stops in coincidence with a complete electron event



Output: electron image





electron image in px

*Projection on x-axis corresponding to
electron energy*



Output: TOF spectra (no coincidences)

Check that peaks in lower subplot are picked correctly!

Eventually play with 'height' in VG_figures.ion_hist line 53



Output: PEPICO matrix (final result)



Output: Files

- *_PEPICO_matrix*.dat text file containing the final PEPICO matrix (columns = TOF, rows = electron energy)
- *_PEPICO_matrix*.png .png image of the PEPICO matrix figure

The 3 numbers at the end of the file name give the TOF upper limit in μs , the binning of the TOF spectrum in ns, and the binning of the electron image energy axis in px.

_9_2_2.dat contains therefore a (45,4500) matrix, for which the energy axis runs from 0 to 90 px in steps of 2 px, and the TOF axis runs from 0 to 9 μs in steps of 2 ns.

Putting together the kinetic energy groups

PEPICO_complete_sp.py

Puts together the PEPICO matrices for the different energy regions to obtain the complete matrix (overlapping regions are averaged)

Input needed in the first section

path to folder containing the input matrices

identifier of the mesurement

parameters used in the previous analysis

work function

path and file to a reference PES

energy calibration factor for PEPICO


central energies of the different regions are retrieved automatically from the Energy .txt file

```
1  # -*- coding: utf-8 -*-
2  """
3  Created on Fri Jan  7 09:32:27 2022
4
5  @author: holzme33
6  """
7  import numpy as np
8  import os
9  import matplotlib.pyplot as plt
10 from matplotlib.gridspec import GridSpec
11 import time
12
13 # PATH = "B:\\analysis\\BSn\\PEPICO\\BSn_168\\"
14 PATH = "B:\\analysis\\BSn\\PEPICO\\BSn_168_fac0.65\\"
15 # PATH = "B:\\analysis\\TBMA\\PEPICO\\TBMA_088_tot\\"
16
17 ide = 'BSn_168'
18 MaxTof = '9'
19 tof_bin = '2'
20 x_step = '2'
21
22 # work function in eV
23 wf = 4.1
24
25 # reference PES (no coincidences)
26 path_ref = "B:\\VG\\rawdata\\"
27 file_ref = "BSn_165Sum"
28 DT_on = 1 # was the reference PES recorded with drift tube on?
29
30 px2eV = 0.85 # calibration factor px to eV
31 DT_shift = 0.6 # energy shift with drift tube on (eV)
32
33 # kinetic energies at which the measurement was done (eV) are retrieved automatically
34 path_en = "B:\\PEPICO\\rawdata\\"
35 file_en = "Energy_+idet+.txt"
36 energy = np.loadtxt(path_en+file_en, skiprows=1, usecols=(2))
37 central_en = np.flip(np.unique(energy))
38
39
40
41 #=====
42 # no need to change anything below this line
43 #=====
44
45 files = []
46 # get the filenames of the raw data saved as .mat
47 for file in os.listdir(PATH):
```


Check parameters, calibration, etc.


Compares the mass-integrated electron spectra for the individual energy regions to the pure photoelectron spectrum.

- Adjust drift tube shift if needed
- Adjust px2eV calibration parameter if needed



Combined PEPICO matrix

- Plots the combined PEPICO matrix and compares the mass-integrated electron spectrum to the reference PES
- Saves the figure as .png
- Saves the matrix and the kinetic energy axis, and the TOF axis in .txt files
- 'cut' (line 85) specifies how many points are dismissed for each energy region at the high energy side
- 'z_max' (line 88) defines the upper limit of the z axis (lower limit is by default 0)



Contact

Address questions, suggestions for extension, bug reports etc. to:

fabian.holzmeier@imec.be